

INCIDENTS | Add Patient



Note: All patient information can be added at any time, before or after patient creation.

01 Apply and Scan Wristband

- Scan wristband to add a patient
- Adding patients by tapping **NEW** increases the risk of a duplicate

02 Set Triage Condition (Color)

03 Add Patient Demographics

- Scan driver's license (DL), or enter manually
- Rapid: **Unknown** name, estimated **Age**, **Gender**

04 Take Photos (Optional)

- Face, injury, DL, belongings, family member, etc.

05 Chief Complaint (Optional)

- Use voice-to-text to describe major injury, vitals, intervention

06 Create Patient

- If you **Select Destination** prior to hitting **Create**, you can immediately alert the receiving facility

