



KNOWLEDGE CHECK

Pulsara MANAGER for EMS & Ambulance

Questions and Answers:

(Correct answers are noted in **BOLD** print)

1. You can access Pulsara MANAGER...
A. from the Pulsara website.
B. from the Pulsara mobile app.
2. The data captured through Pulsara communication can be exported into .CSV and .PDF formats.
A. TRUE
B. FALSE
3. Which of the following statements are TRUE regarding Pulsara passwords.
(Select all that apply.)
A. You can require end users to change their passwords at set intervals.
B. You can require end users to log in using their username and password at set intervals.
C. You can view an end user's password from their user profile.
D. You can reset an end user's password from their user profile.
4. Message templates, such as "CPR in Progress" allow users to efficiently communicate common messages with the simple tap of a button.
A. TRUE
B. FALSE
5. Setting up EMS Units allows end users to choose the unit they are assigned to during their shift.
A. TRUE
B. FALSE

6. Pulsara has integrated with ZOLL to allow ECG transmission from the ZOLL _____ monitor to Pulsara.
- A. ZOLL E-Series
 - B. ZOLL X-Series**
 - C. ZOLL M-Series
7. Pulsara provides Acuity Level selection for which of the following Pulsara modules?
- A. STEMI
 - B. STROKE
 - C. Cardiac Arrest
 - D. Trauma**
 - E. Sepsis
 - F. General Patient**