




# How do I create a general patient case – Apple Devices?


This tutorial shows how an EMS user can start  
and alert the hospital for a general patient

# General Patient



Tap "General" to create a new General case.

# Required Fields



## At Patient Time

Setting the 'At Patient Time' will establish the first medical contact time.

## Patient Name

If patient's name is unknown please enter a value based on current department protocols. An example would be "John Doe".

## Destination Hospital

Choose the destination hospital from the list provided. Only hospitals using Pulsara within the region will be available..


## Required Fields

All fields in red are required to start the case.

## Scene Departure Time

Once scene departure time is set, the distance to the hospital and ETA will be auto calculated using GPS.

# Additional Information



## Optional Patient Information

Age and Gender are optional patient demographics. NOTE: These fields can be important identifiers when the patient arrives to the ED.

## Images

A photo of the Driver's License, ECG, med list, monitor, etc... can be added to the case by tapping the Images/ECG button.

## Chief Complaint

Chief Complaint is a required field to activate a case. Free text all critical information pertaining to the patient so that your hospital partners can prepare the appropriate resources prior to arrival.


# Activating the Case

## Alerting the hospital

Let the hospital know that you are on your way with one simple tap!

Tap 'ALERT' to notify the care team that a new STEMI is in route.

NOTE: Don't forget you can continue to update the case after you have tapped 'Alert'



## Vital Signs

Vital signs are optional. Tapping "+ Vitals" will open the vital signs card and vital information can be entered. These can also be entered after the case is alerted.

## Labs

Typically labs are not drawn in the field. One one that is commonly measured is the blood glucose level. Tapping "+ Labs" will open the field where blood glucose can be entered. These can also be entered after the case is alerted.